



Safe Hygiene at All Levels

OPTUS[®]

DR. NANOXA



+91 8169224938



surgical@rossari.com



ROSSARI BIOTECH LIMITED



Over the last 2.5 decades, Rossari has positioned itself as a leading Specialty Chemicals Manufacturer providing customized solutions to customers' requirements through a diversified product portfolio.

Our Key Segments



Rossari Biotech Limited

Home Personal Care & Performance Chemicals, Textile Specialty Chemicals, Animal Health & Nutrition, Pet Care, BRPL, Performance Additives



Unitop Chemicals Private Limited

(Subsidiary of Rossari Biotech Limited)
Agrochemicals, Oil & Gas, Pharma, Constructions, Paint & Coatings



Tristar Intermediates Private Limited

(Subsidiary of Rossari Biotech Limited)
Personal care, Healthcare Ingredients, Preservatives

Our 4 Pillars of Chemistry

Surfactants

Acrylic Chemistry

Silicones

Enzymes

8 Manufacturing Facilities

435+ Distributors

4220+ Products

Exporting to **33+** Countries & Present all over India

354,100 MTPA Capacity

23+ Industries Served



Ø ZDHC MRSL V3.1 Level 3

(An ISO 9001:2015 & 14001:2015 Certified Company)
201 A - B, 2nd Floor, Akruti Corporate Park, L.B.S Marg, Next to GE Gardens, Kanjurmarg (W), Mumbai - 400078, India

+91 8169224938

info@rossari.com



www.rossari.com

www.unitopchemicals.com

www.tristarintermediates.org

Scan this code to know more about us



HAND SANITIZER

ROSS BACTOKILL

Alcohol Based Effective Handrub

Composition:

Active Ingredients	Quantity
2-Propanol IP	46% v/v
1-Propanol IP	31% v/v
Mecetronium Ethyl Sulphate	0.20% w/w
Other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Application:

- Hand hygiene in hospitals, clinics and other healthcare settings, commercial kitchens, food processing units, educational institutes.
- Preparatory hand disinfection for surgery.

Dilution:

Ready to use.

Activity Time:

- Hygienic hand rub* - 30 secs.
- Surgical hand preparation* - 1.5 to 3 mins.

Certified by:



Available in: 100ml, 500ml, 5ltr

MSDS:



ROSS CIDAL

Isopropyl Alcohol And CHG-Based Handrub

Composition:

Active Ingredients	Quantity
Chlorhexidine Gluconate Solution IP	2.5% v/v (equivalent to Chlorhexidine Gluconate 0.5% w/v)
Isopropyl alcohol IP	70% v/v
Colors and other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Application:

- Hand hygiene in hospitals, clinics and other healthcare settings, commercial kitchens, food processing units, educational institutes.
- Preparatory hand disinfection for surgery.

Dilution:

Ready to use.

Activity Time:

- Hygienic hand rub* - 30 secs.
- Surgical hand preparation* - 1.5 to 3 mins.

Certified by:



Available in: 100ml, 500ml, 5ltr

MSDS:



ROSS CIDAL-E

Ethyl Alcohol And CHG Based Handrub

Composition:

Active Ingredients	Quantity
Chlorhexidine Gluconate Solution IP	2.5% v/v (equivalent to Chlorhexidine Gluconate 0.5% w/v)
Ethyl Alcohol (Ethanol) IP	70% v/v
Other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeastidal.

Application:

- Hand hygiene in hospitals, clinics and other healthcare settings, commercial kitchens, food processing units, educational institutes.
- Preparatory hand disinfection for surgery.

Dilution:

Ready to use.

Activity Time:

- Hygienic hand rub* - 30 secs.
- Surgical hand preparation* - 1.5 to 3 mins.

Certified by:



Available in: 100ml, 500ml, 5ltr

ROSS SANPRO ULTRA

Hand Sanitizer-WHO Recommended Formula

Composition:

Active Ingredients	Quantity
Ethanol IP	80% v/v
Glycerol IP	1.45% v/v
Hydrogen peroxide (H ₂ O ₂) IP	0.125% v/v
Colour Fragrance	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeastidal.

Application:

- Hand hygiene in hospitals, clinics and other healthcare settings, commercial kitchens, food processing units, educational institutes.

Dilution:

Ready to use.

Activity Time:

- Hygienic hand rub* - 30 secs.

Certified by:



Available in: 100ml, 500ml, 5ltr

ANTISEPTICS

ROSSEPT

CHG & Cetrimide Based Antiseptic Cleansing Solution

Composition:

Active Ingredients	Quantity
Strong Cetrimide Solution	0.3% v/v (equivalent to cetrimide IP 3.0% w/v)
Chlorhexidine Gluconate Solution IP	0.30% v/v
Isopropyl Alcohol IP	5% v/v
Benzyl Benzoate IP	0.005% v/v
Purified Water, approved color, perfume	q.s. to 100% v/v

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Applications:

- Cleaning minor burns, wounds, cuts and abrasions.
- First-aid kit.
- General cleaning and disinfection of floors.

Dilution:

- **First-aid** - 15ml of Rossept in 300ml of water.
- **Bathing** - 5ml of Rossept in 10ltr of water.
- **Shaving** - 15ml of Rossept in 500ml of water.

Activity Time:

60 Secs.

Certified by:



Available in: 100ml, 500ml, 5ltr

MSDS:



ROSSEPT DT

Chloroxylenol Based Antiseptic Liquid

Composition:

Active Ingredients	Quantity
Chloroxylenol IP	4.8% w/v
Terpineol BP	9.0% v/v
Alcohol (denatured) eq. to absolute alcohol	13.1% v/v
Other ingredients	Excipients, colours as caramel - q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Applications:

- Cleaning minor burns, wounds, cuts and abrasions.
- First-aid kit.
- General cleaning and disinfection of floors.

Dilution:

- **First-aid** - 15ml of Rossept DT in 300ml of water.
- **Bathing** - 5ml of Rossept DT in 10ltr of water.
- **Shaving** - 12.5ml of Rossept DT in 500ml of water.

Activity Time:

60 Secs.

Certified by:



Available in: 100ml, 500ml, 5ltr

MSDS:



DERMIPREP 5/10

Povidone Iodine Prep-Solution

Composition:

Active Ingredients	Quantity
Povidone-Iodine IP	5 % / 10% w/v (equivalent to 0.5% / 0.10% w/v available iodine)
Other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Application:

- Disinfection of skin surface before any invasive procedure.
- Treatment of infectious wounds, cuts and abrasions.
- Pre and post operative skin disinfection.

Dilution:

Ready to use.

Activity Time:

Immediate after application.

Certified by:



Available in: 100ml, 500ml, 5ltr

DERMIPREP 7.5

Povidone Iodine Surgical Scurb

Composition:

Active Ingredients	Quantity
Povidone-Iodine IP	7.5 % w/v (equivalent to 0.75% w/v available iodine)
Other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Application:

- Treatment and prevention of infections.
- Antiseptic hand wash.
- Pre and post operative surgical hand wash.

Dilution:

Ready to use.

Activity Time:

- Hygienic hand wash* - 30 secs.
- Surgical hand disinfection* - 1.5-3 mins.

Certified by:



Available in: 100ml, 500ml, 5ltr

ROSS CIDAL 4%

Chlorhexidine Gluconate 4% Surgical Hand Scrub

Composition:

Active Ingredients	Quantity
Chlorhexidine Gluconate IP	4% w/v
Colour-Ponceau 4R	q.s

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeastecidal.

Application:

- Surgical hand wash prior to surgery for healthcare professionals.
- Whole body bath (antiseptic wash).
- Routine use.

Dilution:

Ready to Use.

Activity Time:

1.5-3 mins.

Certified by:



Available in: 100ml, 500ml, 5ltr

ALCOHOL SWAB

Composition:

Active Ingredients	Quantity
Saturated isopropyl alcohol	70% v/v

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeastecidal.

Application:

- Skin antiseptic.
- Skin preparation before injection, insertion of IV cannula.

Certified by:



Available in: Pack of 100 Swabs

SURFACE & ENVIRONMENTAL DISINFECTANT

FABSHIELD

Surface Disinfectant

Composition: 100 gm contains,

Active Ingredients	Quantity
Ethanol IP	10 gm
2-Propanol IP	9 gm
1-Propanol IP	6 gm
Other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Application:

- Emergency sanitation of operation theaters between procedures.
- Disinfection of non critical medical devices.

Dilution:

Ready to use.

Activity Time:

30 Secs.

Certified by:



Available in: 250ml, 5ltr

MSDS:



ROSA DC

Surface and Environment Disinfectant

Composition:

Active Ingredients	Quantity
Polyhexamethylene biguanide hydrochloride IP	5 % w/v
Alkyl dimethyl benzylammonium chloride solution IP	15 % w/v
Didecyl dimethyl ammonium chloride solution IP	3% w/v
Other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal, Sporocidal.

Application:

- Cleaning and disinfection of hard, non porous surfaces.
- Aerial disinfection in operation theatre.

Dilution:

- Surface disinfection - 1% (10 ml/Ltr).
- Fogging - 2% (20ml/Ltr).

Activity Time:

- Fumigation time - 20 mins (use of OT after 30 min of fumigation).

Certified by:



Available in: 100ml, 1ltr, 5ltr

MSDS:



ROSSARGEN

Surface and Environmental Disinfectant

Composition:

Active Ingredients	Quantity
Silver accelerated hydrogen peroxide	11% w/w
Silver	150 PPM
Other ingredients	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeastcidal, Sporocidal.

Application:

- Disinfection of surfaces in healthcare settings.
- Aerial disinfection in operation theater.

Dilution:

- **Surface Disinfection (Spray):** Use 10% (10ml ROSSARGEN + 100ml of D.M. water)
- **Aerial Disinfection:** Use 20% (20ml ROSSARGEN + 100ml of D.M. water)

Activity Time:

- Surface disinfection - 30 mins.
- Fogging - 60 mins.

Certified by:



Available in: 100ml, 1ltr, 5ltr

MSDS:



Rosa DC 256

Surface & Environmental Cleaning

Composition:

Active Ingredients	Quantity
Didecyl dimethyl benzyl ammonium chloride	8.70% w/v
n-alkyl dimethyl benzyl ammonium chloride	8.19% w/v
Excipients	83.11% w/v

Efficacy: Bactericidal, Fungicidal, Virucidal, Mycobactericidal, Sporocidal, Mildewstate.

Application:

- Disinfection of surfaces in healthcare settings.
- Aerial disinfection in operation theater.

Dilution:

- **Surface disinfection** - 0.4% (4ml in 1ltr of water).
- **Aerial Disinfection** - 2% (20ml in 1ltr of water).

Activity Time:

20 mins. (Post drying)

Certified by:



Available in: 100ml, 500ml, 5ltr

FABSHIELD EXTRA

Surface & Equipment Disinfectant-Concentrate

Composition: Each 100gm contains,

Active Ingredients	Quantity
(Ethylenedioxy) dimethanol	14.1 g
Glutaraldehyde	5.0 g

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal, Tuberculocidal.

Application:

- Disinfection of surfaces in healthcare settings.
- Aerial disinfection in operation theater.

Dilution:

- **Surface disinfection** - 0.5% (5ml in 1ltr of water).
- **Aerial Disinfection** - 1% (10ml in 1ltr of water).

Activity Time:

5-15 mins.

Certified by:



Available in: 100ml, 500ml, 5ltr

FABSHIELD SPECIAL

Surface & Equipment Disinfectant-Concentrate

Composition: Each 100gm contains,

Active Ingredients	Quantity
1,6 dihydroxy 2,5 dioxahexane (chemically bound formaldehyde)	11.2 gm
Glutaraldehyde	5.0 gm
Benzalkonium chloride IP	5.0 gm
Alkyl Urea Derviates	3.0 gm

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal.

Application:

- Disinfection of surfaces in healthcare settings.
- Aerial disinfection in operation theater.

Dilution:

- **Critical surface disinfection** - 2% (20ml in 1ltr of water).
- **General cleaning and disinfection** - 0.5% (5ml in 1ltr of water).

Activity Time:

30 mins.

Certified by:



Available in: 100ml, 500ml, 5ltr

INSTRUMENT DISINFECTANT

ROSSALDEHYDE OPA

Ortho-Phthalaldehyde High-Level Instrument Disinfectant Solution

Composition:

Active Ingredients	Quantity
Ortho-phthalaldehyde	0.55% w/w
Purified Water IP	q.s.
Other ingredient	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal, Tuberculocidal.

Application:

Disinfection of heat sensitive and semi critical medical devices.

Activity Time:

- Disinfection (Temperature 20°-40°C) - 15-20 mins.
- Sterilization - 10 Hrs.

Certified by:



Available in: 100ml, 1Ltr, 5Ltr with test stripes

MSDS:



ROSSALDEHYDE

Glutaraldehyde Based Instrument Disinfectant Solution

Composition:

Active Ingredients	Quantity
Glutaraldehyde solution IP	2.45% w/v
Purified water IP	q.s.
Excipientis	q.s.

Efficacy: Bactericidal, Fungicidal, Virucidal, Yeasticidal, Tuberculocidal.

Application:

Disinfection of heat sensitive and semi critical medical devices.

Activity Time:

- Disinfection (Temperature 20°-40°C) - 15-20 mins.
- Sterilization - 10 Hrs.

Certified by:



Available in: 1ltr, 5ltr with test stripes

MSDS:



ROSSAZYME

Enzymatic Instrument Cleaner

Composition:

Multi Protease, Cellulose, Amylase, Lipase, Nonionic Surfactant, Chelating Agents, Rust Inhibitors, pH Regulators, Stabilizers.

Efficacy: Removal of lipids, proteins in mucous, saliva and other bodily fluids.

Application:

- Precleaning of critical and non- critical surgical medical devices.
- Remove stains and organic matter such as dried blood, tissues, mucous, etc.

Dilution:

5ml in 1ltr of water.

Activity Time:

- Heavy soiled - Dipping of instruments 10ml for 8min.
- In dried blood contamination dose will 20ml for 20min with leukworm water.

Certified by:



Available in: 1ltr, 5ltr

MSDS:



SURGICAL DISPOSABLE

LATEX EXAMINATION GLOVES

Application:

- Medical & healthcare settings.
- Pharmaceutical manufacturing.

Features:

- Single use
- Ambidextrous

Available in: 20/ 50/ 100 Pcs per pack

Certified by:



Available Size:

Small / Medium / Large

NITRILE EXAMINATION GLOVES

Application:

- Medical & healthcare settings.
- Food processing areas.
- HORECA industries.

Features:

- Single use
- Ambidextrous

Available in: 20/ 50/ 100 Pcs per pack

Certified by:



Available Size:

Small / Medium / Large

SURGICAL GLOVES

Application:

- Surgical procedures.

Features:

- Single use
- Sterilized by ETO

Certified by:



Available sizes:

• 6.0 • 6.5 • 7.0 • 7.5

DISPOSABLE FACE MASK

Features:

- 3 ply
- With and without filter
- Allergy free
- High quality

AVAILABLE IN	UNIT
Face mask with filter	100 pcs
Face mask without filter	100 pcs



BOUFFANT CAP

Features:

- Fluid repellent
- Comfortable and breathable
- Light-weight
- 8 GSM spunbonded fabric

AVAILABLE IN	UNIT
Bouffant Cap	100 pcs



SHOE COVER

Features:

- Provides secure fit
- Prevents slipping and discomfort.

AVAILABLE IN	UNIT
Shoe Cover	50 Pairs



SURGEON CAP

Features:

- Total head coverage
- Elastic provides secure fit
- Comfortable fit

AVAILABLE IN	UNIT
Surgeon Cap	100 pcs



DRESSINGS & ORTHOCARE PRODUCTS

GAUZE SWAB

Application:

Gauze swab are used for wound dressing, padding, protection from blood spills.

Features:

- High volume of retention
- Free from optical whiteness & foreign matter
- Soft & skin friendly

NON-STERILE

SPECIFICATION	UNIT
5 CM X 5 CM X 8 Ply	100 Pieces / Packet
5 CM X 5 CM X 12 Ply	100 Pieces / Packet
7.5 CM X 7.5 CM X 8 Ply	100 Pieces / Packet
7.5 CM X 7.5 CM X 12 Ply	100 Pieces / Packet
10 CM X 10 CM X 8 Ply	100 Pieces / Packet
10 CM X 10 CM X 12 Ply	100 Pieces / Packet



STERILE

SPECIFICATION	UNIT
5 CM X 5 CM X 12 Ply	20 Pieces / Packet
7.5 CM X 7.5 CM X 12 Ply	20 Pieces / Packet
10 CM X 10 CM X 12 Ply	20 Pieces / Packet



GAUZE THAN

Application:

Used by institutions to furnish gauze swabs, gamjee rolls, combined dressings, etc, for their CSSD.

Features:

- Consistent layers
- Even edges
- Highly absorbent



SPECIFICATION	UNIT	PACKAGING
18 mtr X 90 CM Type 13	Single Piece	50 Packets / Bale
18 mtr X 90 CM Type 17	Single Piece	50 Packets / Bale
18 mtr X 100 CM FII - 550gm	Single Piece	50 Packets / Bale

ABD MOP

Application:

For large heavy draining wounds, prevents residue and leakage.

Features:

- With X-ray detectable thread
- Highly absorbent
- Free from foreign matter and optical whiteners
- Sturdy loop attached firmly to one corner

NON-STERILE

SPECIFICATION	UNIT
11 CM X 33 CM (8 Ply, 12 Ply)	10 Pieces / Packet
25 CM X 25 CM (8 Ply, 12 Ply)	10 Pieces / Packet
30 CM X 30 CM (8 Ply, 12 Ply)	10 Pieces / Packet
25 CM X 40 CM (8 Ply, 12 Ply)	10 Pieces / Packet
15 CM X 15 CM (8 Ply, 12 Ply)	10 Pieces / Packet



STERILE

SPECIFICATION	UNIT
25 CM X 25 CM X 8 Ply	1 Piece
30 CM X 30 CM X 8 Ply	1 Piece

GAMJEE ROLL

Application:

Generally used for primary wound dressing.

Features:

- Excellent absorbency
- 100% cotton



SPECIFICATION	UNIT	PACKAGING
10 CM X 3 MTR	10 CM X 12 Roll	144 Packets / Carton
15 CM X 3 MTR	15 CM X 12 Roll	144 Packets / Carton

CAST PADDING

Application:

- Provides cushion between the skin and plaster.



SPECIFICATION	UNIT	PACKAGING
10 CM X 3 MTR	10 CM X 12 Roll	144 Packets / Carton
15 CM X 3 MTR	15 CM X 12 Roll	144 Packets / Carton

ROLLED BANDAGE

Application:

- Used to secure dressing.

SPECIFICATION	PACKAGING
5 CM X 3 MTR X 12 Units	144 Packet / Case
7.5 CM X 3 MTR X 12 Units	96 Packet / Case
10 CM X 3 MTR X 12 Units	72 Packet / Case
15 CM X 3 MTR X 12 Units	48 Packet / Case



TRANSPARENT PU TAPE

Application:

- Dressing and device securement.

Features:

- Hypoallergenic



SPECIFICATION	UNIT	PACKAGING
1.25 CM X 9 MTR / 5 MTR	0.5 inch / 24 Roll	60 Box
2.5 CM X 9 MTR / 5 MTR	1 inch / 12 Roll	60 Box
5 CM X 9 MTR / 5 MTR	2 inch / 6 Roll	60 Box
7.5 CM X 9 MTR / 5 MTR	3 inch / 4 Roll	60 Box

MICROPORE TAPE

Application:

- Dressing and device securement.

Features:

- Hypoallergenic



SPECIFICATION	UNIT	PACKAGING
1.25 CM X 9 MTR	0.5 inch / 24 Roll	60 Box
2.5 CM X 9 MTR	1 inch / 12 Roll	60 Box
5 CM X 9 MTR	2 inch / 6 Roll	60 Box
7.5 CM X 9 MTR	3 inch / 4 Roll	60 Box

ELASTIC ADHESIVE BANDAGE

Application:

Suitable for post operation dressing.

Features:

- Good adhesion to secures heavy dressings
- Porous construction & adhesive coating
- Highly elastic backing



Available in: 10CM X 4/6 MTR

CREPE BANDAGE

Application:

- Help to relieving joint pain & muscle strain, spasm.

Features:

- Better elasticity
- Premium and breathable material



FAST EDGE

SPECIFICATION	PACKAGING
6 CM X 4 MTR	320 Packet / Case
8 CM X 4 MTR	280 Packet / Case
10 CM X 4 MTR	240 Packet / Case
15 CM X 4 MTR	160 Packet / Case

ADULT HYGIENE PRODUCTS

DIAPER

Features:

- High absorbtion
- Spill proof
- Fastenable with tape
- Soft and breathable
- Special weightness indicator



Available Size: M,L,XL

ADULT PULL-UP PANTS

Features:

- High absorbtion
- Spill proof
- Soft and breathable
- Special weightness indicator
- Stretchable waist band



Available Size: M,L,XL & XXL

UNDERPAD

Features:

- Soft
- Absorbent and durable

Application:

- Prevents staining and soiling from bodily fluids and helps keep beds/ linens dry and clean.



Available Size: 90 cm X 60 cm

MEDICAL DEVICES

DIGITAL THERMOMETER

Features:

- Measures actual temperature in °C and ° Farenhite
- Auto switch off

Application:

Measure body temperature.



B.P. MONITOR

Features:

- Fully automatic blood pressure monitor
- Intelligent pressurization
- Arrhythmia detection
- One touch operation

Application:

Measures systolic and diastolic blood pressure and the heart rate pulse.



NEBULIZER

Features:

- Easy medication
- Convenient
- Easy to use
- Fit for all

Applicable Conditions:

Asthma, bronchitis, cold, cough, sore throat, pharyngitis and other respiratory diseases.



AIR MATTRESS

Features:

- Highly effective against bed sores.
- Anti-decubitus
- Air mattress with compressor alternating pressure system.

Application:

Ideal for bed-ridden home patients.





The greatest
wealth is health

Buzil Rossari Private Limited

Unit No F-6, Shree Rajlaxmi Textile Park,
Opp. Delhi Darbar Dhaba, Old Bhiwandi Nashik Road, Village Pogaon, Tal.- Bhiwandi, Dist.- Thane, India.
Email : surgical@rossari.com, Customer Care : **+91 8169224938**